



# Pluripotent Stem Cell Marker Antibodies

# Anti-SSEA-3, SSEA-4, SSEA-1 Monoclonal Antibodies

Stage-Specific Embryonic Antigens (SSEAs) are glycan antigens known to vary in their expression levels across the stages of embryonic development and cell differentiation. SSEA-1, also called Lewis X Antigen, is expressed on the surface of mouse embryonic stem cells (ES cells) and induced pluripotent stem cells (iPS cells) and is often used as a mouse pluripotent stem cell marker. SSEA-3 and SSEA-4 are the Globo-series Glycosphingolipids Gb<sub>5</sub> and SialyIGb<sub>5</sub>, and are localized on the cell surface of human embryonic cancer cells (EC cells), ES cells, iPS cells and Muse (Multilineage-differentiating Stress Enduring) cells, where they are used as markers. SSEA-4 has also been shown to be overexpressed in many cancers, suggesting a correlation with disease progression. (2)

TCI has a wide line-up of mouse monoclonal antibodies, synthetic glycans and glycan conjugates for SSEA-3, SSEA-4, SSEA-1.

### Anti-SSEA-3 (Gb<sub>5</sub>) Monoclonal Antibody

0.1mg/vial [A3729]

Isotype: Mouse IgG3 Specificity: SSEA-3 (Gb<sub>5</sub>)

# Anti-SSEA-4 (SialyIGb<sub>5</sub>) Monoclonal Antibody

0.1mg/vial [A3742]

Isotype: Mouse IgG2b

Specificity: SSEA-4 (SialylGb<sub>5</sub>)


## **Anti-Lewis X Monoclonal Antibody**


0.1mg/vial [A2578]

Isotype: Mouse IgM

Specificity: Lewis X, SSEA-1

# Binding specificity of A3729 and A3742







■ Gb<sub>5</sub> Ceramide [G0592]■ Gb<sub>4</sub> Ceramide [G0580]

Gb<sub>3</sub> Ceramide [G0624]

LactosylceramideEthanol

ELISA plates were coated with the indicated glycolipids. Subsequently, anti-SSEA-3 antibody [A3729] or anti-SSEA-4 antibody [A3742] were added to the wells and allowed to incubate. Primary antibodies were detected using appropriate secondary antibodies.

A3729 and A3742 show high specificity for each antigen.

# **Pluripotent Stem Cell Marker Antibodies:** Anti-SSEA-3, SSEA-4, SSEA-1 Monoclonal Antibodies

Oligosaccharide Structure	Product		
	Glycan	Conjugate	Antibody
SSEA-3 (Gb₅)	Gb <sub>5</sub> Ceramide 1 mg / 5 mg [G0592]  Gb <sub>5</sub> -PrN <sub>3</sub> 1 mg / 5 mg [G0599]  Globopentaose Please contact us. [G0434]	HSA-Gb <sub>5</sub> 0.1mg/vial [H1777] Gb <sub>5</sub> -β- <i>p</i> NP 5mg [G0355]	Anti-SSEA-3 (Gb <sub>5</sub> ) Monoclonal Antibody 0.1mg/vial [A3729]
SSEA-4 (SialylGb <sub>5</sub> )	<b>Gb₅ Ceramide</b> 1mg / 5mg [G0580]		Anti-SSEA-4 (SialylGb <sub>5</sub> ) Monoclonal Antibody 0.1mg/vial [A3742]
Lewis X/SSEA-1	SSEA-1-PrNH <sub>2</sub> Please contact us. [S0946]	HSA-Lewis X 0.1mg/vial [H1719]  Lewis X-Magnetic Beads 0.2mg/vial [L0381]	Anti-Lewis X Monoclonal Antibody 0.1mg/vial [A2578]

References 1) J. K. Henderson, P. W. Andrews, et al., Stem Cells 2002, 20, 329. https://doi.org/10.1634/stemcells.20-4-329

2) K. Sivasubramaniyan, G. Niederfellner, H. J. Bühring, et al., Glycobiology 2015, 25, 902. https://doi.org/10.1093/glycob/cwv032

#### Ordering and **Customer Service**

#### TCI AMERICA

: 800-423-8616 / 503-283-1681 Tel Tel : 888-520-1075 / 503-283-1987 E-mail: Sales-US@TCIchemicals.com

#### TCI EUROPE N.V.

Tel : +32 (0)3 735 07 00 : +32 (0)3 735 07 01 E-mail: Sales-EU@TCIchemicals.com

#### **TCI Deutschland GmbH**

: +49 (0)6196 64053-00 : +49 (0)6196 64053-01 E-mail: Sales-DE@TCIchemicals.com

#### Tokyo Chemical Industry UK Ltd. Tokyo Chemical Industry (India) Pvt. Ltd.

Tel : +44 (0)1865 78 45 60 E-mail : Sales-UK@TCIchemicals.com

#### 梯希爱(上海)化成工业发展有限公司

Tel :800-988-0390 / 021-67121386 : 021-6712-1385 E-mail: Sales-CN@TCIchemicals.com

: 1800 425 7889 / 044-2262 0909 E-mail: Sales-IN@TCIchemicals.com

#### TOKYO CHEMICAL INDUSTRY CO., LTD.

Tel : +81 (0)3-5640-8878

E-mail:globalbusiness@TCIchemicals.com

<sup>•</sup> Chemicals itemized in this brochure are for research and testing use only. Please avoid use other than by chemically knowledgeable professionals. • Information such as listed products and its specifications and so on are subject to change without prior notice. • The contents may not be reproduced or duplicated in whole or in part without permission of Tokyo Chemical Industry Co., Ltd.