


TCIMAIL

No. 197 | AUTUMN 2024



Contents

Research Article 2

Organotellurium-Mediated
Radical Polymerization (TERP)

Shigeru Yamago
Institute for Chemical Research, Kyoto University

Chemistry Chat 15

When a Teacher Is at a Loss in High School

Hiroyuki Onuki
Toyo University Keihoku Senior High School

New Products Information 18

Versatile Trifluoromethylthiolating Reagent: TTST

Diurea and Triurea Derivatives for the Synthesis of Highly Crystalline COFs

Spiroheptadiene Reagents to React with Maleimide to Modify Antibodies

Staining Kit for Staining Viable and Dead Cells

Modulators of the ERK Signaling Pathway

Retinoic Acid Receptors Agonist

Research Article

Organotellurium-Mediated Radical Polymerization (TERP)

Shigeru Yamago

Institute for Chemical Research, Kyoto University

Keywords: controlled radical polymerization, reversible deactivation radical polymerization, tellurium, hyperbranched polymer

1. Introduction

Controlled radical polymerization, also known as living radical polymerization, has made significant advances during the last three decades and has revolutionized the development of polymeric materials.^[1,2] While atom transfer radical polymerization (ATRP) and reversible addition-fragmentation chain transfer polymerization (RAFT) are the most popular methods, the present authors have independently developed organotellurium-mediated radical polymerization (TERP) and have implemented it in society through industry-academic collaborations.^[3-6] TERP fulfills five important technological requirements, which are the keys to its success: 1) high versatility for applicable monomers with different reactivities, 2) high compatibility toward polar functional groups, 3) high flexibility in copolymer synthesis, 4) high diversity in end-

group transformations, and 5) high reliability in largescale synthesis. In this article, the features of TERP related to 1-4) will be discussed. Until recently, controlled radical polymerization was limited to synthesizing polymers with linear structures. However, the authors have recently succeeded in synthesizing hyperbranched polymers with controlled structures by developing monomers that selectively induce branching. This result will also be discussed. The controlled radical polymerization methods reported so far cannot completely prevent the termination reaction. Therefore, it should be noted that IUPAC recommends not using the term “living radical polymerization” for this method but instead using the term “reversible deactivation radical polymerization (RDRP).^[7]

2. Polymerization mechanism and polymerization conditions

The key mechanism of RDRP is reversible radical formation from a radical precursor $P-X$ (P = polymer, X = capping group or atom), called the dormant species (Figure 1a).^[2,3] A P radical generated from the dormant species reacts with a monomer to form a chain-elongated P radical, which is then deactivated to a dormant species. This is different from the formation of dead polymers by the termination reaction in conventional radical polymerization in that the dormant species can regenerate radicals again after its deactivation. In other words, the

propagating end is “living” in the form of the dormant species. The equilibrium is heavily shifted to the dormant species, which reduces the concentration of P radicals. Thus, the probability of the termination reaction occurring is significantly reduced (but not completely eliminated, as noted earlier). Furthermore, if the deactivation occurs faster than or as fast as the propagation reaction, all polymer chains grow at approximately equal rates. The molecular weight of the resulting polymer is predetermined by the monomer/chain transfer agent (or

initiator) ratio with low dispersity (narrow molecular weight distribution), which is a characteristic feature of a living polymerization.

Two mechanisms are known to realize the equilibrium shown in **Figure 1a**: the dissociation-combination (DC) mechanism (**Figure 1b**) and the degenerative chain transfer (DT) mechanism (**Figure 1c**). In the DC mechanism, the formation of P and X radicals by C-X homolysis is the activation step, and the recombination of the formed radical pairs is the deactivation step. In the DT mechanism, on the

other hand, the substitution reaction between a P radical and a dormant species results in simultaneous deactivation of the P radical and activation of the dormant species. For example, ATRP proceeds by the DC mechanism (a high-valent metal formed by the reduction of a dormant species, e.g., Cu(II)X_2 , is the X -radical equivalent), and RAFT proceeds by the DT mechanism. TERP proceeds predominantly by the DT mechanism, but the DC mechanism can also be involved^[8,9] depending on the polymerization conditions for TERP.




Figure 1. a) General mechanism of RDRP, b) DC mechanism, c) DT mechanism, and d) mechanism involving ditelluride. P and P' stand for polymer and X for capping group; in TERP $X = \text{TeR}^1$.

TERP can be carried out under three conditions: 1) a thermal condition by heating the TERP chain transfer agent (CTA) and monomers,^[4,10] 2) a thermal condition by adding an azo initiator to condition 1),^[8] and 3) a photo condition by irradiating a mixture of TERP CTA and monomers with light.^[11,12] In condition 1), the initiating radicals are supplied by thermal dissociation of the C-Te bond of the dormant species, i.e., the DC mechanism, and polymerization proceeds mainly by the DT mechanism. The condition is suitable for polymerizing styrene and dienes, which have a low C-Te bond dissociation energy for the dormant species and a slow propagation rate.^[4,11] Under condition 2), thermal decomposition of the azo initiator provides the initiating radicals, and polymerization

proceeds only by the DT mechanism. This condition is highly versatile since the polymerization temperature can be set according to the decomposition temperature of the azo initiator and the reactivity of the monomer. Under condition 3), the initiating radicals are supplied by photocleavage of the C-Te bond. Since it is important to keep a low concentration of radicals, low-intensity light, such as that from an LED with a power of several watts, is suitable.^[11,13] In addition, polymerization can be carried out at low temperatures, even below room temperature, and the progress of polymerization can be controlled by turning the light on and off (**Figure 2**). Therefore, this condition has unique advantages.




Figure 2. Time control of polymerization by turning the light on and off

Note that the addition of a ditelluride is necessary to achieve high polymerization control in the polymerization of methacrylates.^[10] Ditelluride reacts with the propagating P radical and deactivates it to the dormant species about 100 times faster than the DT mechanism (Figure 1d).^[14] In addition, the generated tellanyl radical ($R^1Te\cdot$) is virtually unreactive toward the monomer

while reacting with the dormant species to regenerate the P radical and the ditelluride. This fast and selective deactivation allows for better control of polymerization in methacrylate polymerization. In contrast, polymerization of acrylate is greatly inhibited by the presence of ditelluride, since activation of the polyacrylate-dormant species by the tellanyl radical is insufficient.

3. TERP chain transfer agent (CTA)

The structures of typical CTAs of TERP are shown in Figure 3. To achieve controlled polymerization, the initiation reaction must occur as fast as or faster than the propagation reaction. Since the radicals generated from CTAs **1**, **2**, and **3** are stable and sufficiently reactive toward the monomer, these CTAs generally give excellent results.^[4,10,15] Among these, **1** having a polymethacrylate-end structure is the most frequently used CTA. On the other hand, **4** with a polyacrylic-end structure^[16] and **5** with a poly(vinyl ether)-end structure^[17] are less efficient than **1** for control because the radicals generated from them are less stable than that from **1**. However, polymerization using these CTAs can be reasonably well controlled.

The substituent effect of R^1 on tellurium is small.^[18] Aromatic groups such as phenyl group show higher reactivity for the DT reaction than alkyl groups such as methyl and butyl groups, and thus, in principle, exhibit better polymerization control. However, since the difference in control is small, there is no need for extra caution in most cases.

As for the ester substituent R^2 , methyl esters and ethyl esters (**1a**, **1b**) are often used in homogeneous polymerization. On the other hand, by condensing the carboxylic acid of **1c** with various amines, various CTAs with functional groups can be synthesized, such as **6** with amino acids and multivalent CTA **7**. This allows for variation in the initiating end structure of the polymer

(Figure 3b).^[19] Furthermore, treatment of **1c** with a base, such as NaOH, generates water-soluble carboxylate **1d** for use in emulsion polymerization in water.^[20] In addition, by immobilizing the CTA using silyl groups, **1e**^[21] can be used for graft polymerization from a surface.

1~5 can be synthesized in a single step by tellanyl

anions and corresponding organohalides. They are sensitive to oxygen but can be stored for a long time under an inert atmosphere. CTAs **1**, **2**, and **3** are now available from TCI, and we hope that this will expand the use of TERP.




Figure 3. Chemical structures of a) typical CTAs and b) functional CTAs for TERP

4. Polymer engineering with TERP

4.1 Synthesis of homopolymers under homogeneous polymerization. The most significant feature of TERP is its high monomer versatility: monomers with different reactivities, from conjugated to non-conjugated monomers, can be polymerized in a controlled manner using the same CTA, resulting in a polymer with a molecular weight (M_n) close to the theoretical value with narrow dispersity (D) (Figure 4).^[4,8,10,11,21–28] This is the advantage of TERP,

in contrast to ATRP and RAFT where the ligand of the catalyst and CTA must be appropriately selected according to the monomer. The monomer conversion usually reaches over 90%, which is also synthetically attractive. Even in the polymerization of ethylene, in which activation of the dormant species is most difficult, polymerization can be controlled under low ethylene pressure and mild thermal conditions.^[27]




Figure 4. Typical examples of polymerization using TERP. The structure of the monomers and the number-average molecular weight (M_n) and dispersity (D) of the resulting polymers are shown. *Ditelluride is added in the polymerization of methacrylate monomers.

TERP has excellent functional group compatibility, including a highly reactive isocyanate group.^[8,11,16,25,29] Bronsted and Lewis acids are also compatible. Therefore, simultaneous control of stereochemistry and molecular weight is possible by TERP of acrylamides in the presence of a lanthanide triflate as a catalyst, which is effective

for stereo-controlled polymerization of acrylamides (Figure 5).^[30–32] The advantages of TERP have also been demonstrated in the synthesis of polymer monoliths and surface-initiated polymerization.^[33–35]




Figure 5. Dual control of stereochemistry and molecular weight

4.2. Emulsion TERP. TERP can be carried out by emulsion polymerization in water due to the high tolerance of organotellurium compounds to water. Microemulsion polymerization with hydrophobic TERP CTAs, such as **1B**,^[36,37] polymerization-induced self-association (PISA) with macro-CTAs having hydrophilic oligomeric chains,^[38–40] and *ab initio* emulsion polymerization using water-soluble CTA **1d** in the presence of a surfactant have been reported (Figure 6).^[20,41] In the latter case,

polymerization proceeds by irradiation with either low-intensity visible light or by thermal conditions. Traditionally, photopolymerization has been considered difficult because light does not penetrate emulsions. We believe that photo-induced emulsion polymerization becomes possible due to the high permeability of visible light and high photosensitivity of the organotellurium compounds.^[19]




Figure 6. a) Polymerization of MMA by emulsion polymerization and b) solution behavior before and after polymerization

*Brij® is a trademark of Croda Americas LLC.

Since the emulsion polymerization condition suppresses the termination reaction due to the segregation effect, high-molecular-weight polymers can be synthesized using emulsion TERP. In particular, the use of **1Dd** and **1Ed** with hydrophilic ethylene glycol groups instead of **1Ad** with a Me group as the R¹ group is effective, and polymerization of butyl acrylate (BA)^[41] and styrene^[42]

gives the corresponding homopolymers with M_n over 500,000 and 1,000,000, respectively, and low dispersity. The high monomer conversion (>90%) is observed even in the synthesis of these ultra-high molecular weight polymers, and this is advantageous in practical use, including the synthesis of block copolymers shown in the next section.

4.3. Synthesis of Copolymers. RDRP can be used to synthesize block copolymers by sequential addition of monomers, similar to other living polymerization methods, but the synthesis usually suffers from the order of monomer addition. However, TERP has fewer limitations than other RDRP methods. For example, when two monomers are selected from the representative conjugated monomers, i.e., styrene, acrylate, and methacrylate, all possible AB-diblock copolymers can be synthesized regardless of the order of monomers used.^[10] The synthesis of block copolymers from

conjugated and non-conjugated monomers is usually highly challenging due to the different reactivities of the monomers. However, even such a combination is possible with TERP, and several block copolymers have already been synthesized,^{[22,43][25]} including a block copolymer consisting of *N*-vinylformamide (NVF) and *N,N*-dimethylacrylamide (**Figure 7a**).^[28] Since the formyl group of PNVF can be selectively hydrolyzed under mild conditions, copolymers with polyvinylamine as a block can be synthesized without affecting the amide groups.^[28]




Figure 7. Typical results of block copolymer synthesis: a) block copolymer composed of conjugated and non-conjugated monomers, b) block copolymer with high molecular weight, and c) block copolymer with conversion of mechanism from radical to cationic polymerization.

High-molecular-weight polymers can be obtained not only by homopolymerization but also by block copolymerization using emulsion polymerization (**Figure 7b**).^[20,41] Since the macro-CTA obtained in the

first polymerization is dispersed in water as particles, there is no need to homogenize the macro-CTA and the second monomer, unlike in homogeneous conditions. Therefore, block polymerization can be started simply by

adding the second monomer after the conversion of the first monomer.

Vinyl ethers and α -olefins do not homopolymerize but copolymerize with acrylate monomers. TERP shows excellent polymerization control in such copolymerizations. For example, alternating copolymers with a controlled structure can be synthesized by copolymerization of vinyl ethers and acrylic monomers by using an excess amount of vinyl ethers (**Figure 7c**).^[44,45] Furthermore, since the


4.4. Conversion of polymer ends. Since organotellurium compounds are excellent precursors not only for carbon-centered radicals but also for carbanions and carbocations, various transformations of the dormant end are possible (**Figure 8a**).^[4] Such transformations usually proceed with high efficiency, thus various functional groups can be introduced at the polymer ends.

The most popular reaction is the radical-mediated reduction. Not only tin hydride, a typical reducing agent,^[4] but also thiols and telluroles, which are generally less efficient reducing agents, can be used due to the high reactivity of the organotellurium compounds.^[48,49] The polymer end can be selectively deuterated by employing the corresponding deuterated reducing agents. The chain reaction proceeds by irradiation with low-intensity light or addition of an azo initiator. Radical-mediated allylation

copolymer has an alkoxy group at the propagating end, the addition of a Lewis acid, such as $\text{BF}_3 \cdot \text{OEt}_2$, initiates living cationic polymerization of the remaining vinyl ether, resulting in the formation of the corresponding block copolymer. Controlled copolymerization of α -olefins and acrylate monomers in the presence of a Brønsted acid also takes place, and the acid increases the feed of α -olefins into the copolymer.^[46,47]

using allylstannane is also possible.^[4]

As the concentration of radicals generated from the dormant species increases with increasing light intensity, radical-radical coupling reactions take place by irradiation with high-intensity light, such as from a high-pressure mercury lamp. For example, using functionalized TEMPO derivatives as the coupling partner, various functional groups can be introduced to the polymer end.^[11,50–52] Diene-mediated dimerization also takes place when photo-activation of the dormant species is carried out in the presence of dienes, even those with functional groups.^[53,54] For example, after coupling using a diene with hydroxy groups, the introduced hydroxy groups can be used for further synthetic transformations, such as the polymerization of lactide, giving a miktoarm polymer with a controlled structure (**Figure 8b**).




Conversion of the dormant end to a carbanion is possible by a transmetalation reaction, even if the dormant polymers possess polar functional groups, such as poly(meth)acrylates and poly(meth)acrylamides, due to the high reactivity of the organotellurium group.^[4,55,56] The resulting carbanions can react with various electrophiles,

including carbonyl compounds and carbon dioxide, to introduce polar functional groups, which are difficult to introduce by radical reactions. Electrophilic activation is also possible, and the transformation using Friedel-Crafts reactions has already been reported.^[57]

5. Controlled synthesis of hyperbranched polymers

Hyperbranched polymers (HBPs) have characteristic properties such as smaller hydrodynamic volume, lower intrinsic viscosity, and a larger number of terminal substituents than linear polymers. Therefore, HBPs have attracted much attention in both basic and applied research.^[58–61] Self-condensing polymerization using an inimer, a term coined for a monomer with an initiator function, has been used to synthesize HBPs under radical

polymerization.^[62–67] However, control of the three-dimensional (3D) structure, including molecular weight, distribution, and branch structure, cannot be achieved except under special conditions. In contrast, the authors have developed a new method by designing vinyl telluride **8** as a branch-inducing monomer for the synthesis of HBP **9** with a controlled 3D structure (Figure 9a).^[49]



[1Ab] ₀ /[8] ₀	3	15	63
generation	2	4	6
M_n (theo)	41,200	39,600	34,600
M_n (SEC)	42,100	25,600	9,800
M_n (MALS)	53,700	57,300	53,900
\bar{D}	1.55	1.71	1.99

Figure 9. a) Synthesis of HBP using vinyl telluride **8** and b) ideal schematic structure of product and polymerization results (using methyl acrylate [MA] as monomer)

The key feature of **8** is that the strength of its C-Te bond changes significantly after it has reacted as a monomer to form **10**. Since generation of the vinyl radical from **8** is unfavorable, **8** does not serve as an initiator. However, after **8** has acted as a monomer, the resulting **10** can generate stable alkyl radicals. Sequential activation of the C-Te bonds from **10** leads to the branched structure in **9** through intermediate **11**. Indeed, theoretical calculations suggest that the C-Te bond dissociation energies for **10** and **11** are about 50 kJ mol⁻¹ lower than that for **8**.

A characteristic feature of this method is that the number of branches and the branching density can be controlled arbitrarily by changing the ratio of CTA (**1**), **8**, and acrylate monomers. For example, controlling the “generation,” which is the hierarchical structure present in dendrimers, can be controlled by the **8/1** ratio, whereas the molecular weight can be controlled by the monomer/**1** ratio. **Figure 9b** shows the ideal schematic structure of HBPs when varying the **8/1** ratio to 3, 15, and 63. The

table in **Figure 9b** also summarizes the polymerization results when setting the monomer/**1** ratio to 500. The molecular weight determined by SEC becomes smaller than the theoretical molecular weight as the number of branches increases, while the absolute molecular weight determined by the multi-angle laser light scattering (MALS) detector is close to the theoretical value. This is consistent with the fact that the hydrodynamic volume of the polymer is reduced due to the branched structure. Note that the dispersity (*D*) of the synthesized HBPs has a larger value than that of linear polymers. Recent simulations have shown that this is due to the distribution of the number of branches.^[68] In addition, the controlled synthesis of hyperbranched polystyrene^[50] and hyperbranched MMA^[69] becomes possible by changing the structure of the vinyl telluride. We are currently working to apply this method to the fabrication of polymer materials, and we hope to introduce the results in the near future.

6. Summary

TERP satisfies the required technical elements to create high-value-added polymer materials at a high level. In addition, the structure-controlled synthesis of HBPs presented in the last section opens up a new possibility for

using HBPs to produce polymeric materials with improved and/or new properties. We hope that TERP will become more popular in the future and contribute to society by being used to create superior functional materials.

Acknowledgments

The research described in this article was conducted at the Graduate School of Engineering, Kyoto University, the Graduate School of Science, Osaka City University, and the Institute for Chemical Research, Kyoto University. The author would like to express his deepest gratitude to the late Professor Junichi Yoshida for his support at the beginning of the research, and to the students, researchers, and staff members who performed the experiments.

The author also thanks Associate Professor Masatoshi Tosaka for his support in completing this manuscript. This work was supported in part by the Japan Science and Technology Agency, the Japan Society for the Promotion of Science, the Asahi Glass Foundation, the Mitsubishi Foundation, the Nagase Science and Technology Foundation, the Toray Science Foundation, and Otsuka Chemical Co., Ltd.

References

- [1] C. Chatgililoglu, A. Studer, Eds.: *Encyclopedia of Radicals in Chemistry, Biology and Materials*, John Wiley & Sons, Ltd, (2012).
- [2] K. Matyjaszewski, M. Möller, Eds: *Polymer Science: A Comprehensive Reference*, ScienceDirect, (2012).
- [3] S. Yamago, *Chem. Rev.* **2009**, *109*, 5051.
- [4] S. Yamago, K. Iida, J.-I. Yoshida, *J. Am. Chem. Soc.* **2002**, *124*, 2874.
- [5] S. Yamago, and Y. Nakamura: in *Polymer Science: A Comprehensive Reference*, 10 Volume Set, Elsevier, 227 (2012).
- [6] S. Yamago, and Y. Nakamura, in *The Chemistry of Organic Selenium and Tellurium Compounds* (R. Zvi, F. L. Joel, M. Ilan, and P. Saul, Eds.) Wiley (2012).
- [7] A. D. Jenkins, R. G. Jones, G. Moad, *Pure Appl. Chem.* **2010**, *82*, 483.
- [8] A. Goto, Y. Kwak, T. Fukuda, S. Yamago, K. Iida, M. Nakajima, J.-I. Yoshida, *J. Am. Chem. Soc.* **2003**, *125*, 8720.
- [9] Y. Kwak, A. Goto, T. Fukuda, Y. Kobayashi, S. Yamago, *Macromolecules* **2006**, *39*, 4621.
- [10] S. Yamago, K. Iida, J.-I. Yoshida, *J. Am. Chem. Soc.* **2002**, *124*, 13666.
- [11] S. Yamago, Y. Ukai, A. Matsumoto, Y. Nakamura, *J. Am. Chem. Soc.* **2009**, *131*, 2100.
- [12] S. Yamago, *Bull. Chem. Soc. Jpn.* **2020**, *93*, 287.
- [13] Y. Nakamura, S. Yamago, *Beilstein J. Org. Chem.* **2013**, *9*, 1607.
- [14] Y. Kwak, M. Tezuka, A. Goto, T. Fukuda, S. Yamago, *Macromolecules* **2007**, *40*, 1881.
- [15] S. Yamago, K. Iida, M. Nakajima, J.-I. Yoshida, *Macromolecules* **2003**, *36*, 3993.
- [16] S. Yamago, Y. Nakamura, Y. Ukai, M. Yu, in *ACS Symposium Series* **2015**, 295.
- [17] E. Mishima, T. Yamada, H. Watanabe, S. Yamago, *Chem. Asian J.* **2011**, *6*, 445.
- [18] E. Kayahara, S. Yamago, Y. Kwak, A. Goto, T. Fukuda, *Macromolecules* **2008**, *41*, 527.
- [19] W. Fan, Y. Nakamura, S. Yamago, *Chem. Eur. J.* **2016**, *22*, 17006.
- [20] W. Fan, M. Tosaka, S. Yamago M. F. Cunningham, *Angew. Chem. Int. Ed.* **2018**, *57*, 962.
- [21] S. Yamago, Y. Yahata, K. Nakanishi, S. Konishi, E. Kayahara, A. Nomura, A. Goto, Y. Tsujii, *Macromolecules* **2013**, *46*, 6777.
- [22] S.-I. Yusa, S. Yamago, M. Sugahara, S. Morikawa, T. Yamamoto, Y. Morishima, *Macromolecules* **2007**, *40*, 5907.
- [23] S. Kumar, M. Changez, C. N. Murthy, S. Yamago, J.-S. Lee, *Macromol. Rapid Commun.* **2011**, *32*, 1576.
- [24] A. Kermagoret, Y. Nakamura, M. Bourguignon, C. Detrembleur, C. Jérôme, S. Yamago, A. Debuigne, *ACS Macro Lett.* **2014**, *3*, 114.
- [25] H. Nakagawa, S. Yamago, K. Ishihara, S.-I. Yusa, *Kobunshi Ronbunshu* **2015**, *72*, 335.
- [26] W. Li, M. Tosaka, S. Yamago, *ACS Symp. Ser.* **2018**, *129*.
- [27] Y. Nakamura, B. Ebeling, A. Wolpers, V. Monteil, F. D'Agosto, S. Yamago, *Angew. Chem. Int. Ed.* **2018**, *57*, 305.
- [28] W. Fan, S. Yamago, *Angew. Chem. Int. Ed.* **2019**, *58*, 7113.
- [29] Y. Nakamura, K. Nakanishi, S. Yamago, Y. Tsujii, K. Takahashi, T. Morinaga, T. Sato, *Macromol. Rapid Commun.* **2014**, *35*, 642.
- [30] B. Park, Y. Imamura, S. Yamago, *Polym. J.* **2021**, *53*, 53515.
- [31] Y. Imamura, T. Fujita, Y. Kobayashi, S. Yamago, *Polym. Chem.* **2020**, *11*, 7042.
- [32] B. Park, M. Tosaka, S. Yamago, *Polym. J.* **2021**, *53*, 53533.
- [33] G. Hasegawa, K. Kanamori, K. Nakanishi, S. Yamago, *Polymer* **2011**, *52*, 4644.
- [34] J. Hasegawa, K. Kanamori, K. Nakanishi, T. Hanada, S. Yamago, *Macromolecules* **2009**, *42*, 1270.
- [35] J. Hasegawa, K. Kanamori, K. Nakanishi, T. Hanada, S. Yamago, *Macromol. Rapid Commun.* **2009**, *30*, 986.
- [36] Y. Sugihara, Y. Kagawa, S. Yamago, M. Okubo, *Macromolecules* **2007**, *40*, 9208.
- [37] Y. Sugihara, S. Yamago, P. B. Zetterlund, *Macromolecules* **2015**, *48*, 4312.
- [38] M. Okubo, Y. Sugihara, Y. Kitayama, Y. Kagawa, H. Minami, *Macromolecules* **2009**, *42*, 1979.
- [39] Y. Kitayama, H. Moribe, K. Kishida, M. Okubo, *Polym. Chem.* **2012**, *3*, 1555.
- [40] Y. Kitayama, M. Okubo, *Polym. Chem.* **2016**, *7*, 2573.
- [41] Y. Jiang, W. Fan, M. Tosaka, M. F. Cunningham, S. Yamago, *Macromolecules* **2021**, *54*, 10691.
- [42] Y. Jiang, W. Fan, M. Tosaka, S. Yamago, *ACS Macro Lett.* **2022**, *11*, 1331.
- [43] S.-I. Yusa, S. Awa, M. Ito, T. Kawase, T. Takada, K. Nakashima, D. Liu, S. Yamago, Y. Morishima, *J. Polym. Sci., Part A: Polym. Chem.* **2011**, *49*, 2761.
- [44] E. Mishima, S. Yamago, *Macromol. Rapid Commun.* **2011**, *32*, 893.
- [45] E. Mishima, S. Yamago, *J. Polym. Sci., Part A: Polym. Chem.* **2012**, *50*, 2254.
- [46] E. Mishima, T. Tamura, S. Yamago, *Macromolecules* **2012**, *45*, 2989.
- [47] E. Mishima, T. Tamura, S. Yamago, *Macromolecules* **2012**, *45*, 8998.
- [48] S. Yamago, A. Matsumoto, *J. Org. Chem.* **2008**, *73*, 7300.
- [49] Y. Lu, T. Nemoto, M. Tosaka, S. Yamago, *Nat. Commun.* **2017**, *8*, 1863.
- [50] Y. Lu, S. Yamago, *Angew. Chem. Int. Ed.* **2019**, *58*, 3952.
- [51] X. Li, T. Kato, Y. Nakamura, S. Yamago, *Bull. Chem. Soc. Jpn.* **2021**, *94*, 996.
- [52] S. Yamago, E. Kayahara, H. Yamada, *React. Funct. Polym.* **2009**, *69*, 416.
- [53] Y. Nakamura, T. Arima, S. Tomita, S. Yamago, *J. Am. Chem. Soc.* **2012**, *134*, 5336.
- [54] Y. Nakamura, T. Arima, S. Yamago, *Macromolecules* **2014**, *47*, 582.
- [55] E. Kayahara, S. Yamago, in *ACS Symposium Series* (Matyjaszewski, K. Sumerlin, B. and Tsarevsky, N. Eds.) American Chemical Society, **2012**, 1100, 99.
- [56] E. Kayahara, H. Yamada, S. Yamago, *Chem. Eur. J.* **2011**, *17*, 1018.
- [57] T. Yamada, E. Mishima, K. Ueki, S. Yamago, *Chem. Lett.* **2008**, *37*, 650.
- [58] Y. Deyue, G. Chao, F. Holger, Eds.: *Hyperbranched Polymers: Synthesis, Properties, and Applications*, Wiley, (2011).
- [59] C. Gao, D. Yan, *Prog. Polym. Sci.* **2004**, *29*, 183.
- [60] B. I. Vöit, A. Lederer, *Chem. Rev.* **2009**, *109*, 5924.
- [61] N. Hadjichristidis, M. Pitsikalis, S. Pispas, H. Iatrou, *Chem.*

- Rev.* **2001**, *101*, 3747.
- [62] J. M. J. Frechet, M. Henmi, I. Gistov, S. Aoshima, M. R. Leduc, B. R. Grubbs, *Science* **1995**, *269*, 1080.
- [63] C. J. Hawker, J. M. J Fr, R. B. Grubbs, J. Daot, *J. Am. Chem. Soc.* **1995**, *117*, 10763.
- [64] K. Matyjaszewski, S. G. Gaynor, A. H. E. Mu, *Macromolecules* **1997**, *30*, 7034.
- [65] B. Liu, A. Kazlauciuonas, J. T. Guthrie, S. Perrier, *Macromolecules* **2005**, *38*, 2131.
- [66] S. G. Gaynor, S. Edelman, K. Matyjaszewski, *Macromolecules* **1996**, *29*, 1079.
- [67] J. A. Alfurhood, P. R. Bachler, B. S. Sumerlin, *Polym. Chem.* **2016**, *7*, 3361.
- [68] M. Tosaka, H. Takeuchi, M. Kibune, T. Tong, N. Zhu, S. Yamago, *Angew. Chem. Int. Ed.* **2023**, *62*, e202305127.
- [69] Y. Lu, S. Yamago, *Macromolecules* **2020**, *53*, 3209.

Author Information



Shigeru Yamago, Professor, Institute for Chemical Research, Kyoto University

[Brief history] 1991 Graduate School of Science and Engineering, Tokyo Institute of Technology (Ph. D.); 1991 Assistant Professor, Department of Chemistry, Faculty of Science, Tokyo Institute of Technology; 1995 Assistant Professor, Graduate School of Engineering, Kyoto University; 1997 Associate Professor, Graduate School of Engineering, Kyoto University; 2003 Professor, Graduate School of Science, Osaka City University; 2006 Professor, Institute for Chemical Research, Kyoto University, to present

[Major awards received] Incentive Award for Synthetic Organic Chemistry (2001), The 27th CSJ Award for Creative Works (2010), The 4th Society of Synthetic Organic Chemistry, Award for Functional Materials (2012), The 44th Ichimura Prizes for Science (2012), The 35th Technical Award of the Adhesion Society of Japan (2013), The Commendation for Science and Technology by the Minister of Education, Culture, Sports, Science and Technology (2015), the 43rd Harushige Inoue Award (2018), The Society of Polymer Science, Japan Award (2019), The 76th CSJ Award (2023), Aggarwal Lecture (Cornell University) (2024).

[Specialized Fields] Synthetic Chemistry, Organic Chemistry, Polymer Chemistry

[Contact] yamago@scl.kyoto-u.ac.jp

Related Products

Diphenyl Ditelluride	1g	5g	D2718
Dimethyl Ditelluride	100mg	1g	D6090
Ethyl 2-Methyl-2-(methyltellanyl)propanoate	100mg	1g	E1508
2-Methyl-2-(methyltellanyl)propanenitrile	100mg	1g	M3520
Methyl(1-phenylethyl)tellane	100mg	1g	M3521
2-Methyl-2-(methyltellanyl)propanoic Acid		100mg	M3728

Chemistry Chat

When a Teacher Is at a Loss in High School

Toyo University Keihoku Senior High School

Hiroyuki Onuki

1. Introduction

In the previous chat,¹ I introduced the difficulty of keeping up with the changes in high school chemistry. In this paper, I reluctantly disclose examples of how I am at a loss regarding what to do when students ask me questions

that I cannot answer appropriately. If I had sufficient knowledge and experience, I could promptly satisfy the students.

2. Which is the best method for measuring molecular weight?

In Japanese high schools, students learn how to measure molecular weight (or formula weight) using the following four methods²:

- Gas equation (Dumas method)
- Elevation of boiling point
- Depression of freezing point
- Osmotic pressure

The appropriate method must be selected according to the compound, and these methods are not always realistic because we have to consume a sample in grams (exceptions apply). Moreover, it is difficult to obtain accurate molecular weights because special techniques are required to avoid technical errors. As a result, the experiment is no longer a measurement of molecular weight but a competition involving measurement skills. We often receive requests from students who want to accurately determine molecular and formula weights using

smaller quantities and simpler methods. However, this is not easy with these methods, thus I am at a loss how to respond to students.

In Japanese high school textbooks, mass spectrometry is not described in detail, only appearing in a topic section.² If Japanese students in general courses encountered a mass spectrometer, it would make a great impact on them because of its simplicity and accuracy. On the other hand, in the International Baccalaureate course, students learn practical mass spectrometry.³ To my great surprise, interpretation of electron-impact mass spectra was given in the Common Test for University Admissions in Japan, 2024.⁴ Even though cognitive thinking is important, I will be at a loss how to answer the Exam coverage when students ask me.

3. Why do so many compound names end in “[n]”?

Students and colleagues often ask this question, but I am at a loss because I do not have any answers. The names of organic compounds that are popular in our daily life often terminate with an [n] pronunciation and have the following suffixes:

- Hydrocarbons : -ane, -ene, -yne (e. g., methane, polyethylene, toluene)
- Nitrogen-containing compounds: -in, -ine (e.g., glycine, adenine, caffeine)
- Antibiotics: -in (e.g., streptomycin, penicillin)

Why are there limited suffixes corresponding to the compound classes? If I told the students that suffixes followed the IUPAC nomenclature rule,⁵ they would never

be convinced of my answer. Relying on my own personal knowledge, I do not have any answers to date. I wish to ask experts about this nomenclature.

4. How do I explain the amount of substance?

The definition of “amount of substance” was changed in 2019,^{6a} and we have to change its explanation accordingly. However, an experienced teacher in our school insisted that a conventional explanation would be more understandable to students than the updated one:

(Previous Explanation)^{6b}

The mole is the amount of substance in a system that contains as many elementary entities as there are atoms in 0.012 kilogram of carbon 12; its symbol is “mol”.

The number of ¹²C atoms in 12 g of carbon 12 was experimentally determined as $6.02... \times 10^{23}$ /mol, which is called the Avogadro constant, N_A . The Avogadro constant and standard of mass (kilogram) were determined independently and contain errors.

(Current Explanation)^{6a}

One mole contains exactly $6.02214076 \times 10^{23}$

elementary entities. This number is the fixed numerical value of the Avogadro constant, N_A , when expressed in mol^{-1} .

In other words, we determined that one mole is a collection of exactly $6.02214076 \times 10^{23}$ particles. Consequently, the mass of one mole of carbon 12 was revised to 11.9999999958 g.

Even after the change, the mass of one mole of ¹²C atoms is the same as before, within the range of significant decimal places used in high school chemistry (usually up to three digits). I introduce only the current definition to my students and omit the historical background to avoid confusion and to save class time. We must proceed to the next difficult topic, a stoichiometry calculation. I usually agonize over whether or not I should tell them the history.

5. Conclusion

Every day I encounter topics that are difficult for a high school teacher to explain, but many topics are probably obvious to experts. From an educational standpoint, it is a good way to entrust questions to students

so they can find answers. However, as a chemist, I wish to find answers by myself. At the very least, when I know a thing, to hold that I know it; and when I do not know a thing, to allow that I do not know it.⁷

Amendment/Acknowledgments

The topic of molecular polarity was introduced in the previous article.¹ One of the readers pointed out that it was incorrect to judge the polarity of a liquid flowing out of a burette by whether the flow bent when a charged rod was brought close to it. The reader stated that the

bending of the water flow was based on a phenomenon called “electrostatic induction,” in which water is attracted to charged substances, and is not attributed to the polarity of the molecules. Details are available in reference 8. I express my deep appreciation for this suggestion.

References and Notes

1. H. Onuki, *TCIMAIL* **2024**, 196, 16.
2. T. Tatsumi, *et al.*, Ministry of Education, Culture, Sports, Science and Technology Approved Textbook “Chemistry”, **2022**, Suken Publishing. (Japanese)
3. S. Owen, *Chemistry for the IB diploma, Second Edition*; Cambridge University Press, 2014; p. 530.
4. National Center for University Entrance Examinations, Common Test for University Admissions in Japan in 2024, “Chemistry” (Japanese)
https://www.dnc.ac.jp/albums/abm.php?d=666&f=abm00004707.pdf&n=2024_or_29_kagaku.pdf (accessed: September 14, 2024).
5. International Union of Pure and Applied Chemistry, “Nomenclature”, (2013)
<https://iupac.org/what-we-do/nomenclature/> (accessed: September 6, 2024).
6. (a) R. Marquardt, J. Meija, Z. Mester, M. Towns, R. Weir, R. Davis, J. Stohner, *Pure Appl. Chem.* **2018**, 90, 175.
 (b) Bureau International des Poids et Mesures, *The International System of Units (SI) 9th edition*, 2019.
7. The Analects of Confucius, Weizeng No. 2, 17.
8. M. Ziaei-Moayyed, E. Goodman, P. Williams, *J. Chem. Educ.* **2000**, 77, 1520.

Author Information



Hiroyuki ONUKI, Ph.D.

He graduated the University of Tokyo in 1989 and received his Ph.D. degree from Graduate School of the University of Tokyo in 1994. He has worked in Nippon Suisan Kaisha, Ltd., RIKEN, Tokyo Chemical Industry Co., Ltd. and Junten Junior and Senior High School. He has concurrently served as an adjunct lecturer in Tokyo University of Agriculture and Technology, Tokyo Denki University, Graduate School of Yokohama City University, Rikkyo University, and Nihon University. In 2020, he was appointed as a science teacher in Toyo University Keihoku Senior High School.

His research interests are organic natural product chemistry, instrumental analyses, and chemical education.

New Products Information

Versatile Trifluoromethylthiolating Reagent: TTST




S-(Trifluoromethyl) Trifluoromethanesulfonylthioate (= TTST) (1) Product Number: **T4211**
1g 5g



Trifluoromethylthio group is known to have remarkably large hydrophobic parameter among fluorine-type substituents.¹⁾ Therefore, various methods and reagents have been developed in order to improve properties of target compounds in the field of drug discovery research. Nevertheless, there are still some problems like stability and versatility.²⁾

TTST (1) is a novel trifluoromethylthiolating reagent developed by Umemoto *et al.*³⁾ 1 is easy-to-handle because 1 is stable liquid in air. Furthermore, 1 shows highly versatile reactivity. By choosing suitable reaction conditions, three types of reactive species (cations, anions and radicals) can be generated. Therefore, it can be used for versatile reactions such as electrophilic trifluoromethylthiolation of electron-rich aromatic rings and reactive methylenes, Cu-mediated nucleophilic trifluoromethylthiolation of aryl iodides and radical trifluoromethylthiolation of alkenes. In addition, because ratio of available group (-SCF₃) in 1 is high, it is advantageous that wastes derived from 1 is less the conventional reagents and that 1 has high atom-efficiency.

New Reagent



Conventional Reagents



References

- 1) C. Hansch, A. Leo, S. H. Unger, K. H. Kim, D. Nikaitani, E. J. Lien, *J. Med. Chem.* **1973**, *16*, 1207.
- 2) X.-H. Xu, K. Matsuzaki, N. Shibata, *Chem. Rev.* **2015**, *115*, 731.
- 3) Y. Yang, S. Miraghaee, R. Pace, T. Umemoto, G. B. Hammond, *Angew. Chem. Int. Ed.* **2023**, *62*, e202306095.

Related Products

Triphenylphosphine	25g	100g	500g	T0519
9-Mesityl-10-methylacridinium Perchlorate (= Mes-Acr ⁺ -Me ClO ₄ ⁻)		1g	5g	M1774
N-Methyl-N-(trifluoromethylthio)aniline		200mg	1g	M2595
N-(Trifluoromethylthio)phthalimide		1g	5g	T3143
Silver(I) Trifluoromethanethiolate		1g	5g	S0977
Copper(I) Trifluoromethanethiolate		1g	5g	C1159

Diurea and Triurea Derivatives for the Synthesis of Highly Crystalline COFs



1,1'-(1,4-Phenylene)diurea (1)

Product Number: **P3216**
1g 5g


1,1'-([1,1'-Biphenyl]-4,4'-diyl)diurea (2)

Product Number: **B6577**
1g 5g

1,1',1''-([1,3,5-Triazine-2,4,6-triyl]tris(benzene-4,1-diyl))triurea (3)

Product Number: **T4192**
1g

Covalent organic frameworks (COFs) are crystalline porous materials synthesized by the covalent linkage of organic linkers, and are anticipated to be used as molecular separation or storage, photocatalyst, and electronics materials.¹⁾ Diurea derivatives such as 1,1'-(1,4-phenylene)diurea (**1**), are employed in the synthesis of COFs composed of flexible urea linkages.²⁾ In addition, Zhu and Cooper *et al.* have reported a method for synthesizing β -ketoenamine-type COFs using urea-type COFs as precursors.³⁾ Initially, urea-linked COFs were synthesized by the condensation of urea derivatives including **1** and 1,1'-([1,1'-biphenyl]-4,4'-diyl)diurea (**2**) with 2,4,6-triformylphloroglucinol. Since crystallinity of COFs is improved via error-correction process during synthesis, the highly reversible urea bonds contribute to improving crystallinity. Subsequently, heating in a solvent containing a small quantity of water led to the hydrolysis of the urea bonds, followed by imine condensation, resulting in β -ketoenamine-type COFs. These "reconstructed" COFs were shown to have higher crystallinity, porosity and photocatalytic performance compared to the corresponding directly synthesized COFs.



References

- 1) C. S. Diercks, O. M. Yaghi, *Science* **2017**, 355, 923.
- 2) C. Zhao, C. S. Diercks, C. Zhu, N. Hanikel, X. Pei, O. M. Yaghi, *J. Am. Chem. Soc.* **2018**, 140, 16438.
- 3) W. Zhang, L. Chen, S. Dai, C. Zhao, C. Ma, L. Wei, M. Zhu, S. Y. Chong, H. Yang, L. Liu, Y. Bai, M. Yu, Y. Xu, X.-W. Zhu, Q. Zhu, S. An, R. S. Sprick, M. A. Little, X. Wu, S. Jiang, Y. Wu, Y.-B. Zhang, H. Tian, W.-H. Zhu, A. I. Cooper, *Nature* **2022**, 604, 72.

Related Product

2,4,6-Triformylphloroglucinol

200mg

1g

T3688



Spiroheptadiene Reagents to React with Maleimide to Modify Antibodies

***N*⁶-[(Spiro[2.4]hepta-4,6-dien-1-ylmethoxy)carbonyl]-L-lysine (1)** Product Number: **H1840**
100mg 500mg

4-Nitrophenyl Spiro[2.4]hepta-4,6-dien-1-ylmethyl Carbonate (2) Product Number: **N1251**
100mg 500mg

*N*⁶-[(Spiro[2.4]hepta-4,6-dien-1-ylmethoxy)carbonyl]-L-lysine (1) enables the site-specific modification of proteins, making it useful for the preparation of homogeneous antibody-drug conjugates (ADCs).¹⁾ The spiroheptadiene (SCp) moiety of 1, introduced into a specific site of the protein by genetic code expansion, reacts with the maleimide moiety at a rate of 1.8-5.4 M⁻¹ s⁻¹ with the Diels-Alder reaction to form a hexacyclic ring.²⁾ ADCs obtained by reacting 1-incorporated antibodies with drug-conjugated maleimides exhibit better stability and efficacy than cysteine-maleimide conjugates both in mouse serum and *in vivo*.

We also offer 4-nitrophenyl spiro[2.4]hepta-4,6-dien-1-ylmethyl carbonate (2) as a reagent to introduce SCp moiety.



References

- 1) R. Dudchak, M. Podolak, S. Holota, O. Szewczyk-Roszczenko, P. Roszczenko, A. Bielawska, R. Lesyk, K. Bielawski, *Bioorg Chem.* **2024**, 143, 106982.
- 2) A. H. St. Amant, F. Huang, J. Lin, D. Lemen, C. Chakiath, S. Mao, C. Fazenbaker, H. Zhong, J. Harper, W. Xu, N. Patel, L. Adams, B. Vijayakrishnan, P. W. Howard, M. Marelli, H. Wu, C. Gao, J. Read de Alaniz, R. J. Christie, *Bioconjugate Chem.* **2019**, 30, 2340.

Related Products

<i>N</i> -(4 <i>E</i>)-TCO-L-lysine	25mg	T4126
L-Azidohomoalanine	100mg	A3550
Maleimido-mono-amide-DOTA	50mg	M3434

Staining Kit for Staining Viable and Dead Cells



Live/Dead Cell Staining Kit [for Cell Staining] (1)

Product Number: **L0465**
1kit

Live/dead cell staining kit (1) is a kit product containing a pre-conditioned Calcein-AM solution and a propidium iodide (PI) solution. Calcein-AM is a cell membrane-permeable compound used to stain live cells and is hydrolyzed by intracellular esterases and emits intense green fluorescence.¹⁾ On the other hand, PI is a fluorescent dye used to stain dead cells. PI cannot penetrate the plasma membrane of living cells, but only enters dead cells or cells in the late stages of apoptosis, where it intercalates with DNA and emits red fluorescence.²⁾ 1 can stain live cells with green fluorescence and dead cells with red fluorescence in approximately 15 minutes. By staining, live and dead cells can be accurately detected.




Figure. K562 cells treated with 0-20% DMSO were stained by 1. As the concentration of DMSO increases, there are fewer live cells (green fluorescence) and more dead cells (red fluorescence).

References

- 1) F. L. Miles, J. E. Lynch, R. A. Sikes, *J. Biol. Methods* **2015**, 2, 1.
DOI: <https://doi.org/10.14440/jbm.2015.73>
- 2) A. M. Rieger, B. E. Hall, L. T. Luong, L. M. Schang, D. R. Barreda, *J. Immunol. Methods* **2010**, 358, 81.

Related Products

ATP-Luciferase Cell Viability Assay Solution	10mL	A3519
Intracellular Reactive Oxygen Species(ROS) Detection Assay Kit	1kit	I1265
MTT Solution [for Cell proliferation assay] (1mL×5)	1set	M3353
Resazurin (Ready-to-use solution) [for Cell proliferation assay]	25mL	R0195


Modulators of the ERK Signaling Pathway



SL 327 [Optimized for Cell Culture] (1)	Product Number: L0470 1mg
PD 98059 [Optimized for Cell Culture] (2)	Product Number: P3159 1mg
PD 184352 [Optimized for Cell Culture] (3)	Product Number: P3157 1mg
Sorafenib [Optimized for Cell Culture] (4)	Product Number: O0672 1mg
GW 5074 [Optimized for Cell Culture] (5)	Product Number: G0649 1mg
FR 180204 [Optimized for Cell Culture] (6)	Product Number: F1381 1mg
Insulin solution (4mg/mL) [Optimized for Cell Culture] (7)	Product Number: I1257 5mL
rhEGF (Human, Recombinant) (100µg/vial) (8)	Product Number: R0262 1vial

Extracellular signal-Regulated Kinase (ERK), a member of the mitogen-activated protein kinase (MAPK) family, plays a crucial role in regulating cell growth, proliferation, and differentiation. Dysregulation of the ERK pathway has been implicated in various diseases, including cancer.¹⁻³⁾

Modulation of the ERK signaling pathway is known to regulate various biological processes which play roles in cellular behavior. The most common commercial inhibitors of the ERK signaling pathway are the various small molecule kinase inhibitors listed above (1-6). Additionally, TCI also offers several common activators of the ERK pathway including insulin (7), and EGF (8). TCI offers the majority of these modulators at a special [Optimized for Cell Culture] grade.



About [Optimized for Cell Culture]

The [Optimized for Cell Culture] series comprises products that have undergone mycoplasma and endotoxin testing, ensuring their safety for use in cell biology research. These products are ideal for various applications, including disease modeling and regenerative medicine research.


References

- 1) R. Ullah, Q. Yin, A. H. Snell, L. Wan, *Semin. Cancer Biol.* **2022**, *85*, 123.
- 2) K. Carel, J. L. Kummer, C. Schubert, W. Leitner, K. A. Heidenreich, B. Draznin, *J. Biol. Chem.* **1996**, *271*, 30625.
- 3) K. Rubio, A. J. Romero-Olmedo, P. Sarvari, G. Swaminathan, V. P. Ranvir, D. G. Rogel-Ayala, J. Cordero, S. Günther, A. Mehta, B. Bassaly, P. Braubach, M. Wygrecka, S. Gattenlöhner, A. Tresch, T. Braun, G. Dobreva, M. N. Rivera, I. Singh, J. Graumann, G. Barreto, *Theranostics* **2023**, *13*, 2384.

Retinoic Acid Receptors Agonist



TTNPB (1)

Product Number: T3842
5mg 25mgTTNPB
(1)

TTNPB (1) is a synthetic stilbene retinoid analog of retinoic acid and acts as a retinoic acid receptor (RAR) agonist (Table 1).¹⁾ 1 also enhances reprogramming efficiency in chemically induced pluripotent stem cells (CiPSCs).²⁾


Table 1. Transcriptional activation assay of 1¹⁾

Retinoic acid receptors (RARs)	EC ₅₀ (nM)
RAR α	21
RAR β	4.0
RAR γ	2.4

References

- 1) R. L. Beard, R. A. S. Chandraratna, D. F. Colon, S. J. Gillett, E. Henry, D. K. Marler, T. Song, L. Denys, M. E. Garst, T. Arefieg, E. Klein, D. W. Gil, L. Wheeler, D. M. Kochhar, P. J. A. Davies, *J. Med. Chem.* **1995**, *38*, 2820.
- 2) P. Hou, Y. Li, X. Zhang, C. Liu, J. Guan, H. Li, T. Zhao, J. Ye, W. Yang, K. Liu, J. Ge, J. Xu, Q. Zhang, Y. Zhao, H. Deng, *Science* **2013**, *341*, 651.

Pamphlets of TCI Products




Transition Metal Catalysts

<https://bit.ly/4e0hWDv>



Macrocycles for Host-Guest Chemistry

<https://bit.ly/4dVFSba>



www.TCIchemicals.com



Ordering and Customer Service

TCI AMERICA

Tel : 800-423-8616 / 503-283-1681
Fax : 888-520-1075 / 503-283-1987
E-mail : Sales-US@TCIchemicals.com

TCI EUROPE N.V.

Tel : +32 (0)3 735 07 00
Fax : +32 (0)3 735 07 01
E-mail : Sales-EU@TCIchemicals.com

TCI Deutschland GmbH

Tel : +49 (0)6196 64053-00
Fax : +49 (0)6196 64053-01
E-mail : Sales-DE@TCIchemicals.com

Tokyo Chemical Industry UK Ltd.

Tel : +44 (0)1865 78 45 60
E-mail : Sales-UK@TCIchemicals.com

梯希爱(上海)化成工业发展有限公司

Tel : 800-988-0390 / 021-67121386
Fax : 021-6712-1385
E-mail : Sales-CN@TCIchemicals.com

Tokyo Chemical Industry (India) Pvt. Ltd.

Tel : 1800 425 7889 / 044-2262 0909
E-mail : Sales-IN@TCIchemicals.com

TOKYO CHEMICAL INDUSTRY CO., LTD.

Tel : +81 (0)3-5640-8878
E-mail : globalbusiness@TCIchemicals.com

• Chemicals itemized in this brochure are for research and testing use only. Please avoid use other than by chemically knowledgeable professionals. • Information such as listed products and its specifications and so on are subject to change without prior notice. • The contents may not be reproduced or duplicated in whole or in part without permission of Tokyo Chemical Industry Co., Ltd.