

New

蛋白酶抑制剂混合物

100×Protease Inhibitor Cocktail (EDTA free)

1 vial [P2949]

100×Protease Inhibitor Cocktail

1 vial [P2976]


优势

- 可溶解在1 mL去离子水中作为100倍储备溶液使用。
- 保护蛋白质样品，以免在多种蛋白酶的作用下发生蛋白质水解。
- **P2949** 不含EDTA，使其与金属螯合物亲和层析兼容。

成分

蛋白酶抑制剂	靶蛋白酶	P2949	P2976
AEBSF	Serine protease	○	○
Aprotinin	Serine protease, Esterase	○	○
EDTA	Metalloprotease		○
E-64	Cysteine protease	○	○
Leupeptin	Cysteine protease, Trypsin-like protease	○	○

蛋白酶抑制试验



和我们竞争对手的产品相比，P2949 和 P2976 显示出对多种蛋白酶更高的抑制作用。

相关产品

AEBSF [for Biochemical Research]	100mg / 1g [A2215]
Benzamidinium Hydrochloride [for Biochemical Research]	5g [B3379]
PMSF [for Biochemical Research]	5g / 25g [B3473]
EDTA 2Na Dihydrate [for Biochemical Research]	5g / 25g [D3789]
EDTA 3Na Hydrate [for Biochemical Research]	5g / 25g [T2599]
EGTA [for Biochemical Research]	5g / 25g [E0805]
E-64d [for Biochemical Research]	5mg / 25mg [E1337]
2-Iodoacetamide [for Biochemical Research]	5g [I0741]
1,10-Phenanthroline Monohydrate [for Biochemical Research]	5g [P1826]
Ubemimex Hydrochloride [for Biochemical Research]	25mg [B4000]