

$^1\text{H-NMR}$



Product No. T2654

2,4,6-Triisopropylphenylboronic Acid
(contains varying amounts of Anhydride)

Solvent : CDCl_3

Internal Standard : $\text{Si}(\text{CH}_3)_4$

Measured Temperature : 20.5°C



^{13}C -NMR



Product No. T2654

2,4,6-Triisopropylphenylboronic Acid
(contains varying amounts of Anhydride)

Solvent : CDCl_3

Internal Standard : CDCl_3

Measured Temperature : 21.2°C



TOKYO CHEMICAL INDUSTRY CO., LTD.